# INSTRUCTIONS FOR USE

Ambu® Reusable PEEP 10 Valve

Ambu® Reusable PEEP 20 Valve

# Ambu

## Content

## Page

|               |    |
|---------------|----|
| English.....  | 4  |
| Česky.....    | 13 |
| Dansk.....    | 22 |
| Deutsch.....  | 32 |
| Ελληνικά..... | 43 |
| Español.....  | 54 |
| Suomi.....    | 65 |
| Français..... | 74 |
| Magyar.....   | 84 |
| Italiano..... | 94 |

## Content

## Page

|                  |     |
|------------------|-----|
| 日本語.....         | 104 |
| Nederlands ..... | 113 |
| Norsk.....       | 123 |
| Polski.....      | 132 |
| Português.....   | 143 |
| Slovenčina.....  | 153 |
| Svenska .....    | 163 |
| Türkçe .....     | 173 |
| 中文.....          | 183 |

## **1. Important information – Read before use**

Read these safety instructions carefully before using the Ambu® Reusable PEEP (Positive End Expiratory Pressure) 10 Valve and Ambu® Reusable PEEP 20 Valve (in this document collectively referred to as Ambu Reusable PEEP Valve unless otherwise specified). The instructions for use may be updated without further notice. Copies of the current version are available upon request. Please be aware that these instructions do not explain or discuss clinical procedures. They describe only the basic operation and precautions related to the operation of the Ambu Reusable PEEP Valve. Before initial use of the Ambu Reusable PEEP Valve, it is essential for operators to have received sufficient training and be familiar with the intended use, warnings, cautions, and indications mentioned in these instructions. There is no warranty on the Ambu Reusable PEEP Valve.

### **1.1. Intended use**

The Ambu Reusable PEEP Valve is suitable for use with resuscitators, portable ventilators, and CPAP systems which have airtight connections and can operate with a positive pressure in the patient connection throughout the entire expiratory phase.

#### **1.1.1. Intended patient population**

The Ambu Reusable PEEP Valve is to be used on all age groups from premature infants to elderly.

#### **1.1.2. Intended use environment**

The Ambu Reusable PEEP Valve is intended to be used in various environments where the need for PEEP might arise. This includes both pre-hospital and hospital settings, including MR system rooms (except inside the MR bore during an ongoing MRI procedure).

### **1.2. Indications for use**

Use of the Ambu Reusable PEEP Valve is indicated when achievement of a pressure above ambient in the patient's lungs is needed in order to improve oxygenation of the blood, and to treat different respiratory diseases and symptoms.

### **1.3. Intended user**

Medical professionals trained in airway management and use of the Ambu Reusable PEEP Valve function such as anesthesiologists, nurses, rescue personnel and emergency personnel.

## 1.4. Contraindications

None known.

## 1.5. Clinical benefits

To improve the functional residual capacity and thus optimise gas exchange, oxygenation and lung compliance.

## 1.6. Warnings and Cautions

### WARNINGS

1. Always visually inspect the product after unpacking, assembly and prior to use as defects and foreign matter can lead to no or reduced ventilation of the patient.
2. Do not reuse the product on another patient without reprocessing due to the risk of cross infection.
3. Only to be used by intended users who are familiar with the content of this manual, as incorrect use might harm the patient.
4. Professionals performing the procedure should assess the choice of size in accordance with the patient's specific condition(s), as incorrect use may harm the patient.
5. Always reprocess the Ambu Reusable PEEP Valve if visible residues are left inside the device in order to avoid the risk of infection or malfunction.
6. Always reprocess the Ambu Reusable PEEP Valve after each use in order to avoid the risk of infection.
7. Do not use the product if contaminated by external sources as this can cause infection.
8. Always use the Ambu Reusable PEEP Valve with a manometer as the manometer ensures correct adjustment of PEEP during ventilation. Too high or low PEEP during ventilation may cause barotrauma and hypoxia, respectively.
9. Do not use the Ambu Reusable PEEP Valve after a maximum of 30 times of reprocessing in order to avoid the risk of infection or malfunction of the device.
10. Always reprocess the Ambu Reusable PEEP Valve minimum once a week when used on the same patient in order to avoid the risk of infection.
11. Do not use the product in patients with obstructive lung disease, unless a medical assessment indicates the necessity, as increased airway pressure can accelerate airway irritation and inflammatory response.

12. Do not use the product in patients with a medical condition causing elevated peak and mean airway pressures, unless a medical assessment indicates the necessity, as increased PEEP can lead to barotrauma.
13. Do not use the product in patients suffering from conditions such as cardiogenic shock, myocardial infarction, left heart failure, hypovolaemia or a combination of these conditions resulting in haemodynamic instability unless a medical assessment indicates the necessity, as high PEEP can lead to reduced cardiac output and systemic perfusion.
14. Do not use a product with decreased performance unless a medical assessment indicates the necessity as decreased performance can lead to hypoxia.
15. Do not use the product in patients with unilateral (one-sided) lung disease, bronchopleural fistulae or similar conditions unless a medical assessment indicates the necessity, as this can lead to hyperinflation of the lung tissue reducing healing of injured site.
16. When using supplemental oxygen, do not allow smoking or use of device near open flame, oil, grease, other flammable chemicals or equipment and tools, which cause sparks, due to the risk of fire and/or explosion.

## **CAUTIONS**

1. Please see packaging for more specific information about the expiration date, as the use of an expired device might lead to decreased performance or malfunction of the product.
2. Always keep components from same device together during reprocessing to avoid reassembly of components with different durability leading to the risk of product failure.
3. Be careful not to unscrew the adjustment cap as it can detach. If this happens the cap can be reattached again.
4. U.S. federal law restricts this device to sale by or on the order of a licensed health care practitioner.

### **1.7. Undesirable Side Effects**

Potential undesirable side effects may include lung injury, barotrauma (e.g., pneumothorax), reduced cardiac output, and decreased systemic perfusion.

## 1.8. General notes

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

## 2. Device description

The Ambu Reusable PEEP Valve is a positive end-expiratory pressure valve for use with breathing equipment conforming to EN ISO 5356-1. The product can be supplied with a connector of either ID Ø 30 mm or OD Ø 22 mm. The Ambu Reusable PEEP Valve is delivered non-sterile and must be reprocessed according to section 3.2.

## 3. Use of Reusable PEEP Valve

### 3.1. Operating the Ambu Reusable PEEP Valve

1. Inspect the Ambu Reusable PEEP Valve to make sure it is free of obstructions.
2. Firmly attach the inlet connector of the Ambu Reusable PEEP Valve to the resuscitator, portable ventilator, or CPAP system in accordance with the Manufacturer's Instructions for Use.
3. Rotate the adjustment cap to the desired PEEP setting
  - PEEP 10: from 0 – 10 cmH<sub>2</sub>O (0 – 1.0 kPa).
  - PEEP 20: from 1.5 – 20 cmH<sub>2</sub>O (0.15 – 2.0 kPa).
4. For correct adjustment of the PEEP, a manometer should be connected to the breathing device for monitoring of the PEEP.

### 3.2. Reprocessing: cleaning, disinfection, sterilization

Follow these reprocessing instructions after each use to reduce the risk of cross-contamination.

#### Disassembly

Before reprocessing, manually disassemble the Ambu Reusable PEEP Valve into three components as shown in illustration 1. Turn the housing counterclockwise while holding the inlet connector. Take out the internal valve assembly. No further disassembly is necessary.

Keep components from same device together during reprocessing to avoid reassembly of components with different durability.

## Recommended reprocessing procedures

For complete reprocessing of the Ambu Reusable PEEP Valve, use one of the procedures listed below.

### Recommended reprocessing procedures (select one)


- Manual cleaning followed by chemical disinfection, or
- Manual cleaning followed by sterilization, or
- Automated cleaning, including a thermal disinfection stage, followed by sterilization, or
- Automated cleaning, including a thermal disinfection stage, followed by chemical disinfection.

Table 1: Recommended reprocessing procedures.

Product testing has shown that the Ambu Reusable PEEP Valve is fully functional after 30 full reprocessing cycles, as listed in Table 1.

It is the user's responsibility to qualify any deviations from the recommended cycles and methods of processing, and to monitor that the recommended number of reprocessing cycles is not exceeded.

### Illustration 1:

## Procedures for reprocessing

### MANUAL CLEANING

1. Rinse the components under running cold utility (tap) water to remove gross soil.
2. Prepare a detergent bath using a cleaning detergent solution, e.g. Neodisher® MediClean Forte or equivalent, for the removal of residues of dried and denatured blood and proteins, using the detergent manufacturer's recommended concentration.
3. Fully immerse the components to keep them submerged in the solution per detergent instruction label. During the soak time thoroughly clean the components with a soft brush until all visible soil is removed.
4. Rinse the components with running utility (tap) water for one minute.
5. Dry the components with a clean lint-free, free cloth and compressed air.

### AUTOMATED CLEANING WITH THERMAL DISINFECTION

1. Rinse the components under running cold utility (tap) water to remove gross soil.
2. Place the components onto a manifold rack or in a wire basket inside the washer.
3. Select the appropriate cycle as listed below:

| Stage                | Recirculation time (minutes) | Temperature                 | Detergent type and concentration   |
|----------------------|------------------------------|-----------------------------|--|
| Pre-wash             | 02:00                        | Cold tap water              | N/A  |
| Wash                 | 01:00                        | 43 °C (110 °F)<br>Tap water | Neodisher® MediClean Forte or an equivalent detergent using manufacturer's recommended concentration |
| Rinse                | 01:00                        | 43 °C (110 °F)<br>Tap water | N/A  |
| Thermal disinfection | 05:00                        | 91 °C (196 °F)              | N/A  |
| Dry time             | 07:00                        | 90 °C (194 °F)              | N/A  |

Table 2: Automated cleaning with thermal disinfection cycles.

## **CHEMICAL DISINFECTION**

1. Equilibrate the bath of Cidex OPA, or an equivalent OPA (ortho-phthalaldehyde), disinfectant at the temperature specified in the OPA disinfectant manufacturer's instructions.
2. Ensure the minimum effective concentration (MEC) of the OPA disinfectant using the OPA test strips specified in the OPA disinfectant manufacturer's instructions.
3. Fully immerse the device in the OPA and ensure all air bubbles are removed from the device surface by agitating the device.
4. Allow the device to soak for the time specified by the OPA disinfectant manufacturer's instructions.
5. Thoroughly rinse the device by fully immersing it in purified water, agitating and allowing it to set for a minimum of 1 minute.
6. Repeat step 5 two more times for a total of 3 rinses using a fresh batch of purified water each time.
7. Dry the device using a sterile lint-free cloth.

## **Sterilization**

Sterilize the product using a gravity steam autoclave running a full cycle at 134 – 135 °C (274 – 275 °F) with an exposure time of 10 minutes and a dry time of 45 minutes. Leave the parts to dry and/or cool completely before reassembling the Ambu Reusable PEEP Valve.

## **Inspection of components**

After reprocessing, carefully inspect the components for damage and residuals or excessive wear and replace if necessary. Some methods may cause discoloration of rubber components without affecting the lifetime. In case of material deterioration, e.g. cracking, the components should be discarded and replaced with a new component.

## **Reassembly**

When mounting the internal valve assembly (see illustration 1), ensure that the ends of the central spindle are correctly situated in the hole of the inlet connector and in the housing with adjustment cap, respectively. Turn the housing clockwise while holding the inlet connector for completion of the assembly.

## **3.3. Disposal**

Used products must be disposed of according to local procedures.

## 4. Technical product specifications

### 4.1. Reusable PEEP Valve specifications

|   |   |
|---|---|
| Connector size  | OD Ø 22 mm or ID Ø 30 mm (EN ISO 5356-1)  |
| Adjustment range  | PEEP 10: 0 – 10 cmH <sub>2</sub> O (0 – 1.0 kPa)<br>PEEP 20: 1.5 – 20 cmH <sub>2</sub> O (0.15 – 2.0 kPa) |
| Dimensions  | Length Ø 65 – 78 mm (depending on adjustment),<br>diameter 45 mm  |
| Weight  | Less than 45 g  |
| Operating temperature limit   | -18 °C to +50 °C (-0.4 °F to +122 °F)*  |
| Storage temperature limit   | -40 °C to +60 °C (-40 °F to +140 °F)*   |
| Recommended long term storage in closed packaging at room temperature,<br>away from sunlight. |   |

\* Tested according to EN ISO 10651-4.

### 4.2. MRI Safety Information


The Ambu Reusable PEEP Valves are classified as MR unsafe. It's recommended to have an alternative MR conditional PEEP valve available.

## 5. Troubleshooting

Problem: It is not possible to obtain correct adjustment of the required PEEP level.  
Possible cause: Leak in the system.  
Correction: Reassemble the product.

Problem: Valve noise.  
Possible cause: Insufficient damping.  
Correction: Disassemble the valve, move the valve disc back and forth in relation to the valve spindle (the inner rod).

## 6. Explanation of symbols used

| Symbol indication   | Description   |
|---|---|
|  | Global Trade Item Number                                    |
|  | Lot Number  |
|  | Date of Manufacture<br>Country of Manufacture               |
|  | Medical Device  |
|  | MR Unsafe   |
| <b>Rx Only</b>  | Prescription use only                                       |
|  | UK Responsible Person                                       |
|  | UK Conformity Assessed                                      |
|  | Importer<br>(For products imported into Great Britain only) |

A full list of symbol explanations can be found on [ambu.com/symbol-explanation](https://www.ambu.com/symbol-explanation).

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Ambu and other trademarks are trademarks of Ambu A/S.
Ambu A/S
Ambu Ltd